

# Medical Flowmeter

## Constant Type Flowmeter

**Feature :**

Use crack-resistant tube made of polycarbonate, pure brass flowmeter body with chrome.  
 Flowmeter in 90 degree/ 180 degree are selectable.  
 The flow output accuracy is within  $\pm 10\%$  which meets medical standard.  
 Each unit has been conducted inlet pressure 50psi standard test before shipment.  
 With GMP, ISO13485 & CE certificates.

**Constant Type Flowmeter:**

- Flow rate remains accurate when inlet pressure is variant.  
 Flow Rate: 0~15 lpm  
 Various types of inlet adapter can be selected.



FA-OB-03



FA-AB-03



FA-OB-27  
 CE mark flowmeter  
 MK IV inlet adapter



FA-OB-91  
 AFNOR inlet adapter



FA-OT-01

## Slim Type & Dial Type Flowmeter

Neonate Use: 0 ~3.5 LPM

**Slim Type Flowmeter:**

- Inlet pressure must be fixed at 50 psi to keep flow rate accuracy.  
 Flow Rate: 0~15 lpm  
 Various types of inlet adapter can be selected.



FA-OS-04



FA-OS-24  
 MK IV inlet adapter



FA-AS-02



FA-OM-01



FA-OS-23



FA-OD-02



FA-TD-02

**Dial Flowmeter:**

- Small and handy size, with user-friendly dial design.  
 Flow Rate: 0~15 lpm  
 0~25 lpm  
 Various types of inlet adapter can be selected.  
 Can connect with humidifier.

## Inlet Adapter



## Oxygen Type --- Model & Specification

Inlet Adapter		without adapter	Ohmeda	MK IV	Puritan-Bennett	Kawasaki	DISS	DIN	Chemetron	AFNOR
<b>Constant Type</b>	Singl	FA-OB-03	FA-OB-13	FA-OB-27	FA-OB-33	FA-OB-53	FA-OB-62	FA-OB-71	FA-OB-81	FA-OB-91
	Twin	FA-OT-01	FA-OT-11	FA-OT-21	FA-OT-31	FA-OT-51	FA-OT-62	FA-OT-71	FA-OT-81	FA-OT-91
<b>Slim Type</b>	Singl	FA-OS-04	FA-OS-12	FA-OS-23	FA-OS-32	FA-OS-52	FA-OS-62	FA-OS-71	FA-OS-81	FA-OS-91
	Twin	FA-OM-01	FA-OM-11	FA-OM-21	FA-OM-31	FA-OM-51	FA-OM-62	FA-OM-71	FA-OM-81	FA-OM-91
<b>0-3.5 Slim Type</b>	Singl	FA-OS-05	FA-OS-15	FA-OS-24	FA-OS-33	FA-OS-53	FA-OS-63	FA-OS-72	FA-OS-82	FA-OS-92
<b>Dial Type 15lpm</b>	Singl	FA-OD-02	FA-OD-12	FA-OD-22	FA-OD-32	FA-OD-52	FA-OD-64	FA-OD-72	FA-OD-82	FA-OD-92
	Twin	FA-TD-02	FA-TD-12	FA-TD-22	FA-TD-32	FA-TD-52	FA-TD-64	FA-TD-72	FA-TD-82	FA-TD-92
<b>Dial Type 25lpm</b>	Singl	FA-OD-01	FA-OD-11	FA-OD-21	FA-OD-31	FA-OD-51	FA-OD-63	FA-OD-71	FA-OD-81	FA-OD-91
	Twin	FA-TD-01	FA-TD-11	FA-TD-21	FA-TD-31	FA-TD-51	FA-TD-63	FA-TD-71	FA-TD-81	FA-TD-91

## Air Type --- Model & Specification

Inlet Adapter		without adapter	Ohmeda	MK IV	Puritan-Bennett	Kawasaki	DISS	DIN	Chemetron	AFNOR
<b>Constant Type</b>	Singl	FA-AB-03	FA-AB-13	FA-AB-27	FA-AB-33	FA-AB-53	FA-AB-62	FA-AB-71	FA-AB-81	FA-AB-91
	Twin	FA-AT-01	FA-AT-11	FA-AT-21	FA-AT-31	FA-AT-51	FA-AT-62	FA-AT-71	FA-AT-81	FA-AT-91
<b>Slim Type</b>	Singl	FA-AS-02	FA-AS-12	FA-AS-22	FA-AS-32	FA-AS-52	FA-AS-62	FA-AS-71	FA-AS-81	FA-AS-91
	Twin	FA-AM-01	FA-AM-11	FA-AM-21	FA-AM-31	FA-AM-51	FA-AM-62	FA-AM-71	FA-AM-81	FA-AM-91
<b>0-3.5 Slim Type</b>	Singl	FA-AS-03	FA-AS-13	FA-AS-23	FA-AS-33	FA-AS-53	FA-AS-63	FA-AS-72	FA-AS-82	FA-AS-92