

Oxygen Regulator with Flowmeter - Single Stage Style

Precision control to reach accurate consumption of gas.

VSY Series



VSW Series



VSC Series



- Piston type with fixed pressure output
- Diaphragm type with outside adjustable pressure output
- Without flowmeter design for central Oxygen Supply
- Output pressure 50psi
- 2" Pressure Gauge
- Crack-resistant polycarbonate tubes in transparent display
- Flow output accuracy within $\pm 10\%$ fulfilled medical standard
- Regulator's body and Flowmeter's body are all Brass chromed plated in for long time use and less maintenance issue.
- CGA 870 (Pin Index) / CGA 540 (Bull- Nose) / CGA 992 (BS) inlet connection can be selected

Inlet



VSW-103



VSW-305

Inlet: CGA 870

Inlet: CGA 540



VSY-203



VSW-218

Inlet: CGA 992

Inlet: CGA 992 with hand-tight wheel



ISO 13485:2003

GMP