


ACARE Click Style Oxygen Regulator -

Model: VST-AM1, VST-AM2



VST series are advanced oxygen regulators which can regulate accurate flow output. The lightweight Aluminum body and all-brass pressure chamber can provide maximum safety and high performance. The recessed adjoining body are well integrated with easy turned flow knob which is clear view in flow precise adjustment.



• Mini size - VST-AM1



Body Length: 9.3 cm

• Mini size - VST-AM2



Body Length: 9.3 cm

Model	Flow	Inlet Connection
Mini size		
VST-301	0- 8 LPM	CGA 870
VST-AM1	0-15 LPM	CGA 870
VST-AM2	0-15 LPM	CGA 870
VST-313	0-15 LPM	CGA 870
VST-303	0-25 LPM	CGA 870
VST-107	0- 8 LPM	CGA 540
VST-108	0-15 LPM	CGA 540
VST-109	0-25 LPM	CGA 540
Regular size		
VST-304	0- 8 LPM	CGA 870
VST-305	0-15 LPM	CGA 870
VST-311	0-15 LPM	CGA 870
VST-306	0-25 LPM	CGA 870
VST-104	0- 8 LPM	CGA 540
VST-105	0-15 LPM	CGA 540
VST-106	0-25 LPM	CGA 540
Long size		
VST-307	0- 8 LPM	CGA 870
VST-308	0-15 LPM	CGA 870
VST-309	0-25 LPM	CGA 870
2 DISS Model - Regular size		
VST-314	0-15 LPM	CGA 870
VST-312	0-25 LPM	CGA 870

FDA ISO 13485:2003 GMP

ACARE Acare Technology Co., Ltd.

6F.-6, No.5, Wuquan 1st Rd., Xinzhuang Dist., New Taipei City 242, Taiwan

Tel: +886-2-2298-8170

Fax: +886-2-2298-8560

Email: info@acaretech.com

http://www.acaretech.com

Acare-ORC02-20110201