



Reusable SPO2 Probe (AccuSensor Series)

Adult SpO2 Probe Sensor



To collocate with patient monitor and pulse oximeter, Acare offers a wide range of very durable and reusable probe sensors for continuous, non-invasive measurement of arterial oxygen saturation. Each sensor has been conducted strict quality inspection and ensures its accuracy. Every sensor is easy to apply and designed with user-friendly. To be compatible with different brand technology, Acare provides versatile probe sensors for selection. Adult probe is comfort and secure fit with minimal pressure. It's suitable for the Operating Room, Recovery Room, Intensive Care Unit and Wards.



ASANR-D1/ Nellcor Compatible

- Plug type: DB9 pin
- Cable length: 1.0 m
- Latex Free
- Compatible: Nellcor, Colin, Datascope (S&W), Drager, Fukuda Dens
- Finger clip SpO2 probe sensor for Adult
- Available weight upon 40 Kgs



- Features:
 - Comfortable fit for continuous monitoring
 - Accurate and dependable oxygen saturation reading
 - Easy for clean and it's durable design
 - Self-aligning grip and secure use
 - Latex-free
 - 12-Month guarantee

Selection availability-

- BCI compatible : ASABI-D1

Plug type: DB9 pin



- Philips compatible : ASAPS-M1

Plug type: 8 pin



- Nonin compatible : ASANN-D1

Plug type: DB9 pin



- Novamatrix compatible : ASANX-M1

Plug type: 6 pin



- Ohmeda compatible : ASAOA-M1

Plug type: 7 pin



Acare Technology Co., Ltd.

6F.-6, No.5, Wuquan 1st Rd., Xinzhuang Dist., New Taipei City 242, Taiwan

Tel: +886-2-2298-8170

Fax: +886-2-2298-8560

Email: info@acaretech.com

www.acaretech.com

Acare-SA01-20110209