




Reusable SP02 Probe (AccuSensor Series)

Pediatric SpO2 Probe Sensor



To collocate with patient monitor and pulse oximeter, Acare offers a wide range of very durable and reusable probe sensors for continuous, non-invasive measurement of arterial oxygen saturation. Each sensor has been conducted strict quality inspection and ensures its accuracy. Every sensor is easy to apply and designed with user-friendly. To be compatible with different brand technology, Acare provides versatile probe sensors for selection. Pediatric probe is an easy finger-clip and comfortable sensor for most of children with less resistance. It's suitable for the Operating Room, Recovery Room, Intensive Care Unit and Wards.



ASPNR-D1 / Nellcor Compatible

- Plug type: DB9 pin
- Cable length: 1.0 M
- Compatible: Nellcor, Colin, Datascope (S&W), Drager, Fukuda Dens
- Finger clip SpO2 probes sensor for Pediatrics
- Available weight upon 10-40 Kgs



- Features:
 - Keeps accurate result while the child is in minor motion
 - For short-term monitoring and quick spot-checking
 - Easy for clean and it's durable design
 - Self-aligning grip and comfortable
 - Latex-free
 - 12-Month guarantee

Selection availability-

- BCI compatible : ASPBI-D1

Plug type: DB9 pin



- Ohmeda compatible: ASPOA-M1

Plug type: 7 pin



Acare Technology Co., Ltd.

6F.-6, No.5, Wuquan 1st Rd., Xinzhuang City, Taipei County 242, Taiwan

Tel: +886-2-2298-8170

Fax: +886-2-2298-8560

Email: info@acaretech.com

<http://www.acaretech.com>