

Reusable SPO2 Probe (AccuSensor Series)

Silicone Soft-tip Probe Sensor

To collocate with patient monitor and pulse oximeter, Acare offers a wide range of very durable and reusable probe sensors for continuous, non-invasive measurement of arterial oxygen saturation. Each sensor has been conducted strict quality inspection and ensures its accuracy. Every sensor is easy to apply and designed with user-friendly. To be compatible with different brand technology, Acare provides versatile probe sensors for selection. Silicone soft-tip probe is good at life expectancy and the risk of cross-contamination is low in any situation. It's suitable for the Operating Room, Recovery Room, Intensive Care Unit and Wards.



ASSNR-D1 / Nellcor Compatible

- Plug type: DB9 pin
- Cable length: 1.0 M
- Compatible: Nellcor, Colin, Datascope (S&W), Drager, Fukuda Dens
- Silicon soft-tip SpO2 sensor for Adults
- Available weight upon 30Kg



- Features:
 - Accurate and dependable SpO2 reading
 - Quick and Easy spot-checking
 - Easy for clean and it's durable design
 - Secure and comfortable in fitting and fixing
 - Latex-free

Selection availability-

- BCI compatible : ASSBI-D1

Plug type: DB9 pin



- Philips compatible : ASSPS-D1

Plug type: 8 pin



- Nonin compatible : ASSNN-D1

Plug type: DB9 pin

